

Nasopharyngoskop mit LED-Lightstick

Nasopharyngoscope with LED-Lightstick




Arbeitslänge	<i>Working Length</i>	300 mm	300 mm
Durchmesser	<i>Diameter</i>	2,8 mm	3,4 mm
Bildwinkel	<i>Direction of View</i>	70°	70°
Abwinklung	<i>Bending distal end</i>	up/down 150°	up/down 150°
Artikel Nr.	<i>Product No.</i>	50 28 30 NP	50 34 30NP

FLEXIBLE *ENDOSKOPE* 300/3.4 + 300/2.8

- einlegbar
- brilliantes Bild
(4 Beleuchtungslichtleiter)
- korrosionsbeständig
- ergonomisches Design
- 18.000 Pixel


Die Symbiose aus Medizinpraxis
und Ingenieurwissen



FLEXIBLE *ENDOSCOPES* 300/3.4 + 300/2.8

- Insertable
- Crystal clear image (4 light guides)
- Non-corrosive
- Ergonomic design
- 18.000 pixel

The symbiosis of medical
knowledge and engineering
know-how



LED - Lightstick

LED 914065
LED Lightstick
Länge: 6,5 cm
Batterienersatz für
200 Std. Brenndauer

LED 914140
LED Lightstick
Länge: 14 cm
Batterienersatz für
250 Std. Brenndauer



LED - Lightstick

LED 914065
LED Lightstick
Length: 6,5 cm
Battery set for
200 operating hours

LED 914140
LED Lightstick
Length: 14 cm
battery set for
250 operating hours